


# ARCTIC JUMBO SHIPPER

## VAPOR SHIPPER

For the safe transport of biological samples in a vapor-stage.



## FEATURES

- 15 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid



SPECIFICATIONS	JUMBO SHIPPER
LN <sub>2</sub> Capacity (L)	15
Static Hold Time (Days)	21
Evap. Rate (Liters/Day)	0.7
DIMENSIONS	
Necktube (in)	8.55
Overall Dia. (in)	18.5
Height (in)	23
WEIGHT	
Empty (Lbs)	40.5
Charged Weight (Lbs)	67
CAPACITY	
.5 cc Straws on Canes	N/A
.5 cc Straws – Bulk	N/A
1.2/2 ml Vials on Canes – Bulk	966
1.2/2.0 ml Vials in Boxes	500
Blood Bags (4r9953)	10
Number of Canisters	N/A



800.886.2796 | Fax 317. 297.7988  
sales@ICBiomedical.com  
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268  
1031 Cass-White Rd. | Cartersville, GA 30121