


ARCTIC SHIPPER X

VAPOR SHIPPER

For the safe transport of biological samples in a vapor-stage.



FEATURES

- 10 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	ARCTIC SHIPPER X
LN ₂ Capacity (L)	10
Static Hold Time (Days)	14
Evap. Rate (Liters/Day)	0.7
DIMENSIONS	
Necktube (in)	8.55
Overall Dia. (in)	15
Height (in)	23
WEIGHT	
Empty (Lbs)	28.9
Charged Weight (Lbs)	44.9
CAPACITY	
.5 cc Straws on Canes	N/A
.5 cc Straws – Bulk	N/A
1.2/2 ml Vials on Canes – Bulk	966
1.2/2.0 ml Vials in Boxes	500
Blood Bags (4r9953)	10
Number of Canisters	N/A