



Vapor Shippers

ARCTIC 2VS

VAPOR SHIPPER



For the safe transport of biological samples in a vapor-stage.



FEATURES

- 2.4 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	ARCTIC 2VS
LN ₂ Capacity (L)	2.4
Static Hold Time (Days)	17
Evap. Rate (Liters/Day)	0.14
DIMENSIONS	
Necktube (in)	2.25
Overall Dia. (in)	8.75
Height (in)	16.375
WEIGHT	
Empty (Lbs)	8
Charged Weight (Lbs)	12.4
CAPACITY	
.5 cc Straws on Canes	55
.5 cc Straws – Bulk	165
1.2/2 ml Vials on Canes – Bulk	N/A
Number of Canisters	1



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

ARCTIC 4VX

VAPOR SHIPPER



For the safe transport of biological samples in a vapor-stage.



FEATURES

- 3.8 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	ARCTIC 4VX
LN ₂ Capacity (L)	3.8
Static Hold Time (Days)	18
Evap. Rate (Liters/Day)	0.2
DIMENSIONS	
Necktube (in)	2.75
Overall Dia. (in)	8.75
Height (in)	18.5
WEIGHT	
Empty (Lbs)	11
Charged Weight (Lbs)	19
CAPACITY	
.5 cc Straws on Canes	280
.5 cc Straws – Bulk	440
1.2/2 ml Vials on Canes – Bulk	95
Number of Canisters	1



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

ARCTIC 7VS

VAPOR SHIPPER



For the safe transport of biological samples in a vapor-stage.



FEATURES

- 6 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	ARCTIC 7V
LN2 Capacity (L)	6
Static Hold Time (Days)	26
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in)	3.75
Overall Dia. (in)	11.25
Height (in)	20.5
WEIGHT	
Empty (Lbs)	15
Charged Weight (Lbs)	25
CAPACITY	
.5 cc Straws on Canes	390
.5 cc Straws – Bulk	969
1.2/2 ml Vials on Canes – Bulk	110
Number of Canisters	1



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

ARCTIC 20VS

VAPOR SHIPPER



For the safe transport of biological samples in a vapor-stage.



FEATURES

- 6.3 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	ARCTIC 20VS
LN2 Capacity (L)	6.3
Static Hold Time (Days)	28
Evap. Rate (Liters/Day)	0.22
DIMENSIONS	
Necktube (in)	3.75
Overall Dia. (in)	15
Height (in)	26.9
WEIGHT	
Empty (Lbs)	29.5
Charged Weight (Lbs)	40
CAPACITY	
.5 cc Straws on Canes	2730
.5 cc Straws – Bulk	6789
1.2/2 ml Vials on Canes – Bulk	770
Number of Canisters	6



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

ARCTIC JUMBO SHIPPER

VAPOR SHIPPER

For the safe transport of biological samples in a vapor-stage.

FEATURES

- 15 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	JUMBO SHIPPER
LN ₂ Capacity (L)	15
Static Hold Time (Days)	21
Evap. Rate (Liters/Day)	0.7
DIMENSIONS	
Necktube (in)	8.55
Overall Dia. (in)	18.5
Height (in)	23
WEIGHT	
Empty (Lbs)	40.5
Charged Weight (Lbs)	67
CAPACITY	
.5 cc Straws on Canes	N/A
.5 cc Straws – Bulk	N/A
1.2/2 ml Vials on Canes – Bulk	966
1.2/2.0 ml Vials in Boxes	500
Blood Bags (4r9953)	10
Number of Canisters	N/A

ARCTIC SHIPPER X

VAPOR SHIPPER

For the safe transport of biological samples in a vapor-stage.

FEATURES

- 10 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	ARCTIC SHIPPER X
LN ₂ Capacity (L)	10
Static Hold Time (Days)	14
Evap. Rate (Liters/Day)	0.7
DIMENSIONS	
Necktube (in)	8.55
Overall Dia. (in)	15
Height (in)	23
WEIGHT	
Empty (Lbs)	28.9
Charged Weight (Lbs)	44.9
CAPACITY	
.5 cc Straws on Canes	N/A
.5 cc Straws – Bulk	N/A
1.2/2 ml Vials on Canes – Bulk	966
1.2/2.0 ml Vials in Boxes	500
Blood Bags (4r9953)	10
Number of Canisters	N/A