

# LABS80K-AI

## BULK SEMEN CRYOSTORAGE



Designed for the Animal Husbandry Market



High efficiency – very low LN<sub>2</sub> consumption.

- Storage capacity of 622,000 0.5cc straws
- LN<sub>2</sub> storage in either liquid or vapor phase

Single pivot point aluminum turntable design – smooth rotation and ergonomic use.

- Eliminates mechanical component maintenance
- Full height, segmented and labeled turntable provides easy inventory storage access and management
- Convenient access to measure LN<sub>2</sub> level
- Available with optional Liquid Level Controller and double-step

Hinged lockable lid and full-size stainless steel tabletop. Easy operation with a convenient working surface.

LABS80K-AI	SPECIFICATIONS
Liquid Nitrogen Capacity (L)	1269
Number of Canisters	311
Canister Dimensions (in)	2.75 Dia. x 25.0 Length
0.5cc Straws (10/cane)	N/A
0.5cc Straws- Bulk (5 Levels) in Goblets (65mm goblet - 400 straws/goblet)	622,000
0.25cc Straws (28/cane)	N/A
0.25cc Straws - Bulk (5 Levels) in Goblets (65mm goblets - 860 straws/goblet)	1,337,300



800.886.2796 | Fax 317. 297.7988  
 sales@ICBiomedical.com • www.ICBiomedical.com