

TW 38HC

LN₂ STORAGE SMALL FREEZER



Designed for storing large quantities of AI/IVF samples at extremely low temperatures



FEATURES

- Wide neck openings and large capacity provides easy access to stored samples
- Provides sample security and simplified inventory control
- Smaller footprint and lighter full shipping weight
- Numbered index rings and tapered internal canister location for ease of use
- Rugged construction - ribbed high-strength aluminum body, magenformed necktube design and durable paint
- 5-Year vacuum warranty

SPECIFICATIONS	38HC
LN ₂ Capacity (L)	38
Static Hold Time (Days) Liquid	120
Working Time Days (Liquid)	75
Evap. Rate (Liters/Day)	0.31
DIMENSIONS	
Canister Dimensions (In)	2.75 x 11
Neck Diameter (in)	5.0
Overall Height (in)	26.8
Overall Diameter (in)	18.8
WEIGHT	
Weight Empty (Lbs)	38.2
Weight Full (Lbs)	105.9
CAPACITY	
Number of Canisters	10
# of 1.2 & 2.0 ml Vials (5/Cane)	950
# of 1.2 & 2.0 ml Vials (6/Cane)	1140
# of .05 cc Straws (10/Cane)	3200
# of .25 cc Straws (28/Cane)	8960
# of .05 cc Straws Bulk (2 Levels) in Goblet	8000
# of .25 cc Straws Bulk (2 Levels) in Goblet	17,200



800.886.2796 | Fax 317. 297.7988
 sales@ICBiomedical.com • www.ICBiomedical.com

© 2021, IC Biomedical, LLC