

# TW 38VHC

## LN<sub>2</sub> HIGH-CAPACITY SMALL FREEZER

For long-term storage of large quantities at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid



SPECIFICATIONS	38VHC
LN <sub>2</sub> Capacity (L)	38
Static Hold Time (Days)	120
Evap. Rate (Liters/Day)	0.31
DIMENSIONS	
Necktube (in)	5.00
Overall Dia. (in)	18.8
Height (in)	26.8
WEIGHT	
Empty (Lbs)	38.2
Full Weight (Lbs)	104.7
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	3528
.5 cc Straws – Bulk (1 Level)	5760
Number of 1.2 & 2.0 ml Vials	1260



800.886.2796 | Fax 317. 297.7988  
[sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)  
[www.ICBiomedical.com](http://www.ICBiomedical.com)

© 2021, IC Biomedical, LLC