

TW CXR 100

VAPOR SHIPPER

For the safe transport of biological samples in a vapor phase.



FEATURES

- 3.7 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	CXR100
LN2 Capacity (L)	3.7
Static Hold Time (Days)	16
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in)	3.58
Overall Dia. (in)	9.2
Height (in)	19.4
WEIGHT	
Empty (Lbs)	11.7
Charged Weight (Lbs)	18.2
CAPACITY	
No. of 1.2 ml & 2.0 ml vials (5/cane)	85
No. of 1.2 ml & 2.0 ml vials (6/cane)	102
Number of .25 cc Straws (28/cane)	784
No. of .25 cc Straws, Bulk (2 levels) in Goblets	1640
Number of 1/2 straws (10/cane)	280
No. of .5 cc Straws, Bulk (2 levels) in Goblets	800
Number of Canisters	1



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268
1031 Cass-White Rd. | Cartersville, GA 30121

© 2020, IC Biomedical, LLC