

TW CXR500

VAPOR SHIPPER



For the safe transport of biological samples in a vapor phase.



FEATURES

- 6.4 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	CXR500
LN ₂ Capacity (L)	6.4
Static Hold Time (Days)	11
Evap. Rate (Liters/Day)	0.64
DIMENSIONS	
Necktube (in)	8.5
Overall Dia. (in)	15.5
Height (in)	26.9
WEIGHT	
Empty (Lbs)	30
Charged Weight (Lbs)	41.5
CAPACITY	
No. of 1.2 ml & 2.0 ml vials (5/cane)	500
No. of 1.2 ml & 2.0 ml vials (6/cane)	500
Number of .25 cc Straws (28/cane)	1736
No. of .25 cc Straws, Bulk (2 levels) in Goblets	7160
Number of 1/2 straws (10/cane)	1240
No. of .5 cc Straws, Bulk (2 levels) in Goblets	3200
Number of Canisters	N/A



800.886.2796 | Fax 317. 297.7988
 sales@ICBiomedical.com
 www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268
 1031 Cass-White Rd. | Cartersville, GA 30121