

TW 10LD

LN₂ DEWAR

For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	10LD
LN ₂ Capacity (L)	10
Static Hold Time (Days)	45
Evap. Rate (Liters/Day)	0.22
DIMENSIONS	
Necktube (in)	2.0
Overall Dia. (in)	11.4
Height (in)	23.5
WEIGHT	
Empty (Lbs)	14.5
Full (Lbs)	32.3

FEATURES

- 10 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268
1031 Cass-White Rd. | Cartersville, GA 30121

© 2021, IC Biomedical, LLC