

TW 20HC

LN₂ HIGH-CAPACITY SMALL FREEZER



For long-term storage of large quantities at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	20HC
LN ₂ Capacity (L)	20
Static Hold Time (Days)	87
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in)	3.6
Overall Dia. (in)	15.6
Height (in)	24.25
WEIGHT	
Empty (Lbs)	26.4
Full Weight (Lbs)	62
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	1850
.5 cc Straws – Bulk (1 Level)	2940
Number of 1.2 & 2.0 ml Vials	570



800.886.2796 | Fax 317. 297.7988
 sales@ICBiomedical.com
 www.ICBiomedical.com