

TW 35HC

LN₂ HIGH-CAPACITY SMALL FREEZER

For long-term storage of large quantities at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid



SPECIFICATIONS	35HC
LN ₂ Capacity (L)	35
Static Hold Time (Days)	130
Evap. Rate (Liters/Day)	0.27
DIMENSIONS	
Necktube (in)	4.7
Overall Dia. (in)	18.8
Height (in)	26.8
WEIGHT	
Empty (Lbs)	37.9
Full Weight (Lbs)	100.2
CAPACITY	
Number of Canisters	10
.5 cc Straws on Canes	2800
.5 cc Straws – Bulk (1 Level)	4900
Number of 1.2 & 2.0 ml Vials	850



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com