

TW 35LD

LN₂ DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	35LD
LN ₂ Capacity (L)	35
Static Hold Time (Days)	152
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in)	2.5
Overall Dia. (in)	18.8
Height (in)	26.3
WEIGHT	
Empty (Lbs)	35.1
Full (Lbs)	97.4

FEATURES

- 35 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268
1031 Cass-White Rd. | Cartersville, GA 30121