

# TW 50LD

## LN<sub>2</sub> DEWAR

For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	50LD
LN <sub>2</sub> Capacity (L)	50
Static Hold Time (Days)	122
Evap. Rate (Liters/Day)	0.41
DIMENSIONS	
Necktube (in)	2.5
Overall Dia. (in)	18.8
Height (in)	32.4
WEIGHT	
Empty (Lbs)	38.7
Full (Lbs)	127.7

## FEATURES

- 50 Liter LN<sub>2</sub> capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



800.886.2796 | Fax 317. 297.7988  
sales@ICBiomedical.com  
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268  
1031 Cass-White Rd. | Cartersville, GA 30121