

# TW 5LD

## LN<sub>2</sub> DEWAR

For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	5LD
LN <sub>2</sub> Capacity (L)	5
Static Hold Time (Days)	6
Evap. Rate (Liters/Day)	0.77
DIMENSIONS	
Necktube (in)	5.6
Overall Dia. (in)	7.6
Height (in)	17.5
WEIGHT	
Empty (Lbs)	6.9
Full (Lbs)	15.8

## FEATURES

- 5 Liter LN<sub>2</sub> capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



800.886.2796 | Fax 317. 297.7988  
sales@ICBiomedical.com  
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268  
1031 Cass-White Rd. | Cartersville, GA 30121

© 2021, IC Biomedical, LLC