

TW 8XTL

LN₂ EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	8XTL
LN ₂ Capacity (L)	8
Static Hold Time (Days)	80
Evap. Rate (Liters/Day)	0.1
DIMENSIONS	
Necktube (in)	2.00
Overall Dia. (in)	15.6
Height (in)	19
WEIGHT	
Empty (Lbs)	19.6
Full Weight (Lbs)	33.8
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	N/A
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	N/A



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

© 2021, IC Biomedical, LLC