

# TW LS4800

## VIAL STORAGE



For liquid- or vapor-stage vial storage in convenient box-type racks.



SPECIFICATIONS	LS4800
LN <sub>2</sub> Capacity (L)	130
Static Hold Time (Days)	153
Evap. Rate (Liters/Day)	0.85
DIMENSIONS	
Necktube (in)	8.5
Overall Dia. (in)	26.9
Height (in)	35.1
WEIGHT	
Empty (Lbs)	90
Full (Lbs)	312.4

## FEATURES

- Convenient Storage Systems – rack index location ring and internal location fixture, computer compatible box storage for simple inventory management
  - Number of Racks: 6
  - Number of Shelves / Rack: 8
  - 2ml Vial Capacity: 4800
- In vapor-phase storage, temperature at the top of racks is -190°C.
- High strength aluminum body and magneformed necktube design
- Superior vacuum performance with super insulation provides maximum holding times
- Security – Accessory low-level alarm and CS100 controller are available with remote connection capability
- Lockable lids



800.886.2796 | Fax 317. 297.7988  
sales@ICBiomedical.com  
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268  
1031 Cass-White Rd. | Cartersville, GA 30121