

ARCTIC CLASSIC-25

LN₂ DEWAR

For safe storage and dispensing of small accounts of liquid nitrogen.



SPECIFICATIONS	ARCTIC CLASSIC-25
LN ₂ Capacity (L)	24.1
Static Hold Time (Days)	125
Evap. Rate (Liters/Day)	0.20
DIMENSIONS	
Necktube (in)	2.0
Overall Dia. (in)	15.5
Height (in)	22.9
WEIGHT	
Empty (Lbs)	16.5
Full (Lbs)	60.0

FEATURES

- 24.1 Liter LN₂ capacity
- Rugged, lightweight aluminum design
- Designed for long-lasting performance
- 5-year vacuum warranty