

ARCTIC 10

LN₂ SMALL FREEZER



Designed for the safe transport and storage of biological samples.



FEATURES

- Increased capacity
- Rugged lightweight aluminum construction
- Convenient specimen loading/unloading
- Durable exterior finish
- Enhanced necktube strength
- Greater efficiency / low LN₂ loss
- Lockable lid / insert
- 5-Year vacuum warranty

SPECIFICATIONS	ARCTIC 10
LN ₂ Capacity (L)	10
Static Hold Time (Days)	100
Evap. Rate (Liters/Day)	0.1
DIMENSIONS	
Necktube (in)	2.25
Overall Dia. (in)	11.75
Height (in)	23
WEIGHT	
Empty (Lbs)	12
Full Weight (Lbs)	33
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	660
.5 cc Straws – Bulk (1 Level)	879
Number of 1.2 & 2.0 ml Vials	150



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268
1031 Cass-White Rd. | Cartersville, GA 30121