

ARCTIC 20VS

VAPOR SHIPPER



For the safe transport of biological samples in a vapor-stage.



FEATURES

- 6.3 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	ARCTIC 20VS
LN ₂ Capacity (L)	6.3
Static Hold Time (Days)	28
Evap. Rate (Liters/Day)	0.22
DIMENSIONS	
Necktube (in)	3.75
Overall Dia. (in)	15
Height (in)	26.9
WEIGHT	
Empty (Lbs)	29.5
Charged Weight (Lbs)	40
CAPACITY	
.5 cc Straws on Canes	2730
.5 cc Straws – Bulk	6789
1.2/2 ml Vials on Canes – Bulk	770
Number of Canisters	6



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268
1031 Cass-White Rd. | Cartersville, GA 30121