

ARCTIC 20X

LN₂ SMALL FREEZER



Designed for the safe transport and storage of biological samples.



FEATURES

- Increased capacity
- Rugged lightweight aluminum construction
- Convenient specimen loading/unloading
- Durable exterior finish
- Enhanced necktube strength
- Greater efficiency / low LN₂ loss
- Lockable lid / insert
- 5-Year vacuum warranty

SPECIFICATIONS	ARCTIC 20X
LN ₂ Capacity (L)	20
Static Hold Time (Days)	105
Evap. Rate (Liters/Day)	0.19
DIMENSIONS	
Necktube (in)	3.75
Overall Dia. (in)	15
Height (in)	25.5
WEIGHT	
Empty (Lbs)	24
Full Weight (Lbs)	60
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	2100
.5 cc Straws – Bulk (1 Level)	3300
Number of 1.2 & 2.0 ml Vials	630



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268
1031 Cass-White Rd. | Cartersville, GA 30121