

ARCTIC 2VS

VAPOR SHIPPER



For the safe transport of biological samples in a vapor-stage.



FEATURES

- 2.4 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	ARCTIC 2VS
LN ₂ Capacity (L)	2.4
Static Hold Time (Days)	17
Evap. Rate (Liters/Day)	0.14
DIMENSIONS	
Necktube (in)	2.25
Overall Dia. (in)	8.75
Height (in)	16.375
WEIGHT	
Empty (Lbs)	8
Charged Weight (Lbs)	12.4
CAPACITY	
.5 cc Straws on Canes	55
.5 cc Straws – Bulk	165
1.2/2 ml Vials on Canes – Bulk	N/A
Number of Canisters	1



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268
1031 Cass-White Rd. | Cartersville, GA 30121