

# ARCTIC 38-10

## LN<sub>2</sub> SMALL FREEZER



Designed for the safe transport and storage of biological samples.



### FEATURES

- Increased capacity
- Rugged lightweight aluminum construction
- Convenient specimen loading/unloading
- Durable exterior finish
- Enhanced necktube strength
- Greater efficiency / low LN<sub>2</sub> loss
- Lockable lid / insert
- 5-Year vacuum warranty

SPECIFICATIONS	ARCTIC 38-10
LN <sub>2</sub> Capacity (L)	38
Static Hold Time (Days)	135
Evap. Rate (Liters/Day)	0.28
DIMENSIONS	
Necktube (in)	5.00
Overall Dia. (in)	18.5
Height (in)	27.5
WEIGHT	
Empty (Lbs)	31
Full Weight (Lbs)	102
CAPACITY	
Number of Canisters	10
.5 cc Straws on Canes	3200
.5 cc Straws – Bulk (1 Level)	5000
Number of 1.2 & 2.0 ml Vials	1320



800.886.2796 | Fax 317. 297.7988  
sales@ICBiomedical.com  
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268  
1031 Cass-White Rd. | Cartersville, GA 30121