

ARCTIC 38-6

LN₂ SMALL FREEZER



Designed for the safe transport and storage of biological samples.



FEATURES

- Increased capacity
- Rugged lightweight aluminum construction
- Convenient specimen loading/unloading
- Durable exterior finish
- Enhanced necktube strength
- Greater efficiency / low LN₂ loss
- Lockable lid / insert
- 5-Year vacuum warranty

SPECIFICATIONS	ARCTIC 38-6
LN ₂ Capacity (L)	38
Static Hold Time (Days)	135
Evap. Rate (Liters/Day)	0.28
DIMENSIONS	
Necktube (in)	5.00
Overall Dia. (in)	18.5
Height (in)	27.5
WEIGHT	
Empty (Lbs)	31
Full Weight (Lbs)	102
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	3900+
.5 cc Straws – Bulk (1 Level)	6200
Number of 1.2 & 2.0 ml Vials	1200+



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268
1031 Cass-White Rd. | Cartersville, GA 30121