

# ARCTIC 50X

## LN<sub>2</sub> SMALL FREEZER



Designed for the safe transport and storage of biological samples.



### FEATURES

- Increased capacity
- Rugged lightweight aluminum construction
- Convenient specimen loading/unloading
- Durable exterior finish
- Enhanced necktube strength
- Greater efficiency / low LN<sub>2</sub> loss
- Lockable lid / insert
- 5-Year vacuum warranty

SPECIFICATIONS	ARCTIC 50X
LN <sub>2</sub> Capacity (L)	50
Static Hold Time (Days)	125
Evap. Rate (Liters/Day)	0.40
DIMENSIONS	
Necktube (in)	3.75
Overall Dia. (in)	18.5
Height (in)	30.75
WEIGHT	
Empty (Lbs)	38
Full Weight (Lbs)	128
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	2100
.5 cc Straws – Bulk (1 Level)	3300
Number of 1.2 & 2.0 ml Vials	630



800.886.2796 | Fax 317. 297.7988  
sales@ICBiomedical.com  
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268  
1031 Cass-White Rd. | Cartersville, GA 30121