

ARCTIC 7VS

VAPOR SHIPPER



For the safe transport of biological samples in a vapor-stage.



FEATURES

- 6 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	ARCTIC 7V
LN ₂ Capacity (L)	6
Static Hold Time (Days)	26
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in)	3.75
Overall Dia. (in)	11.25
Height (in)	20.5
WEIGHT	
Empty (Lbs)	15
Charged Weight (Lbs)	25
CAPACITY	
.5 cc Straws on Canes	390
.5 cc Straws – Bulk	969
1.2/2 ml Vials on Canes – Bulk	110
Number of Canisters	1



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com

4040 Championship Dr. | Indianapolis, IN 46268
1031 Cass-White Rd. | Cartersville, GA 30121