

# D2000C

## VIAL STORAGE

For liquid or vapor phase vial storage in convenient box-type racks.



SPECIFICATIONS	D2000C
LN <sub>2</sub> Capacity (L)	60
Static Hold Time (Days)	83
Evap. Rate (Liters/Day)	0.72
DIMENSIONS	
Necktube (in)	8.5
Overall Dia. (in)	22
Height (in)	28
WEIGHT	
Empty (Lbs)	55
Charged Weight (Lbs)	162

## FEATURES

- Greater liquid nitrogen efficiency
- Rugged, lightweight aluminum construction
- Large neck opening
- Easy access and retrieval of samples
- Maximum storage capacity
- Durable inventory racks
- Lockable lid