

D4000

VIAL STORAGE

For liquid or vapor phase vial storage in convenient box-type racks.



SPECIFICATIONS	D4000
LN2 Capacity (L)	120
Static Hold Time (Days)	150
Evap. Rate (Liters/Day)	0.8
DIMENSIONS	
Necktube (in)	8.5
Overall Dia. (in)	22
Height (in)	39
WEIGHT	
Empty (Lbs)	71
Charged Weight (Lbs)	285

FEATURES

- Greater liquid nitrogen efficiency
- Rugged, lightweight aluminum construction
- Large neck opening
- Easy access and retrieval of samples
- Maximum storage capacity
- Durable inventory racks
- Lockable lid