

D6000

VIAL STORAGE

For liquid or vapor phase vial storage in convenient box-type racks.



SPECIFICATIONS	D6000
LN ₂ Capacity (L)	165
Static Hold Time (Days)	194
Evap. Rate (Liters/Day)	0.84
DIMENSIONS	
Necktube (in)	8.5
Overall Dia. (in)	26.9
Height (in)	39
WEIGHT	
Empty (Lbs)	121
Charged Weight (Lbs)	410

FEATURES

- Greater liquid nitrogen efficiency
- Rugged, lightweight aluminum construction
- Large neck opening
- Easy access and retrieval of samples
- Maximum storage capacity
- Durable inventory racks
- Lockable lid