

# ARCTIC IC-3D - IC-100D



## LN<sub>2</sub> DEWAR FAMILY

For safe storage and dispensing of small amounts of liquid nitrogen.



## FEATURES

- 3.6 up to 100 Liter capacity
- Rugged, lightweight aluminum design
- Designed for long-lasting performance
- 5-year vacuum warranty

SPECIFICATIONS	IC-3D	IC-5D	IC-6D	IC-10D	IC-20D	IC-35D	IC-35DX	IC-50D	IC-100D
LN <sub>2</sub> Capacity (L)	3.6	5	6.1	10	20.1	35	35	50	100
Static Hold Time (Days)	25	25	30	62	100	159	140	125	250
Evap. Rate (Liters/Day)	.14	.20	.20	.16	.20	.22	.25	.40	.40
DIMENSIONS									
Necktube Dia. (in)	2.25	2.25	2.25	2.25	2.25	2.25	3.75	3.75	3.75
Overall Dia. (in)	8.75	8.75	8.75	11.25	15.00	18.50	18.50	18.50	18.50
Height (in)	16.25	15.75	18.10	23.00	25.50	26.00	26.00	30.75	49.10
WEIGHT									
Empty (Lbs)	6	6	7	12	23	31	31	38	81
Full (Lbs)	13	15	18	33	60	94	94	128	259



4040 Championship Dr. | Indianapolis, IN 46268  
1031 Cass-White Rd. | Cartersville, GA 30121

800.886.2796 | Fax 317. 297.7988

sales@ICBiomedical.com

www.ICBiomedical.com