



LN₂ Dewars

ARCTIC IC-3D - IC-100D



LN₂ DEWAR FAMILY

For safe storage and dispensing of small amounts of liquid nitrogen.



FEATURES

- 3.6 up to 100 Liter capacity
- Rugged, lightweight aluminum design
- Designed for long-lasting performance
- 5-year vacuum warranty

SPECIFICATIONS	IC-3D	IC-5D	IC-6D	IC-10D	IC-20D	IC-35D	IC-35DX	IC-50D	IC-100D
LN ₂ Capacity (L)	3.6	5	6.1	10	20.1	35	35	50	100
Static Hold Time (Days)	25	25	30	62	100	159	140	125	250
Evap. Rate (Liters/Day)	.14	.20	.20	.16	.20	.22	.25	.40	.40
DIMENSIONS									
Necktube Dia. (in)	2.25	2.25	2.25	2.25	2.25	2.25	3.75	3.75	3.75
Overall Dia. (in)	8.75	8.75	8.75	11.25	15.00	18.50	18.50	18.50	18.50
Height (in)	16.25	15.75	18.10	23.00	25.50	26.00	26.00	30.75	49.10
WEIGHT									
Empty (Lbs)	6	6	7	12	23	31	31	38	81
Full (Lbs)	13	15	18	33	60	94	94	128	259



855.750.8191 | sales@ICBiomedical.com

www.ICBiomedical.com

Made In USA

© 2021, IC Biomedical, LLC

ARCTIC CLASSIC-25

LN₂ DEWAR

For safe storage and dispensing of small accounts of liquid nitrogen.



SPECIFICATIONS	ARCTIC CLASSIC-25
LN ₂ Capacity (L)	24.1
Static Hold Time (Days)	125
Evap. Rate (Liters/Day)	0.20
DIMENSIONS	
Necktube (in)	2.0
Overall Dia. (in)	15.5
Height (in)	22.9
WEIGHT	
Empty (Lbs)	16.5
Full (Lbs)	60.0

FEATURES

- 24.1 Liter LN₂ capacity
- Rugged, lightweight aluminum design
- Designed for long-lasting performance
- 5-year vacuum warranty