

# TW 20HC

## LN<sub>2</sub> HIGH-CAPACITY SMALL FREEZER



For long-term storage of large quantities at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	20HC
LN <sub>2</sub> Capacity (L)	20
Static Hold Time (Days)	87
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in)	3.6
Overall Dia. (in)	15.6
Height (in)	24.25
WEIGHT	
Empty (Lbs)	26.4
Full Weight (Lbs)	62
CAPACITY	
Number of Canisters	6
Canister Dimensions (in)	2.75 x 11
1.2ml & 2.0 ml vials (5/cane)	570
1.2ml & 2.0 ml vials (6/cane)	684
.5 cc Straws (10/cane)	1920
.5 cc Straws- Bulk (2 Levels)	4800
.25 cc Straws (28/cane)	5376
.25 cc Straws- Bulk (2 Levels) in Goblets	10740



855.750.8191 | [sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)  
[www.ICBiomedical.com](http://www.ICBiomedical.com)  
 Made In USA