

# TW 35HC

## LN<sub>2</sub> HIGH-CAPACITY SMALL FREEZER

For long-term storage of large quantities at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid



| SPECIFICATIONS                            | 35HC      |
|---|-----------|
| LN <sub>2</sub> Capacity (L)              | 35        |
| Static Hold Time (Days)                   | 130       |
| Evap. Rate (Liters/Day)                   | 0.27      |
| DIMENSIONS                                |           |
| Necktube (in)                             | 4.7       |
| Overall Dia. (in)                         | 18.8      |
| Height (in)                               | 26.8      |
| WEIGHT                                    |           |
| Empty (Lbs)                               | 39.0      |
| Full Weight (Lbs)                         | 101.3     |
| CAPACITY                                  |           |
| Number of Canisters                       | 10        |
| Canister Dimensions (in)                  | 2.64 x 11 |
| 1.2ml & 2.0 ml vials (5/cane)             | 850       |
| 1.2ml & 2.0 ml vials (6/cane)             | 1020      |
| .5 cc Straws (10/cane)                    | 2800      |
| .5 cc Straws- Bulk (2 Levels)             | 7300      |
| .25 cc Straws (28/cane)                   | 7280      |
| .25 cc Straws- Bulk (2 Levels) in Goblets | 16400     |



855.750.8191 | [sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)  
[www.ICBiomedical.com](http://www.ICBiomedical.com)  
Made In USA