

TW 38VHC
**LN₂ HIGH-CAPACITY
SMALL FREEZER**

For long-term storage of large quantities at cryogenic temperatures.


FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	38VHC
LN ₂ Capacity (L)	38
Static Hold Time (Days)	120
Evap. Rate (Liters/Day)	0.31
DIMENSIONS	
Necktube (in)	5.00
Overall Dia. (in)	18.8
Height (in)	26.8
WEIGHT	
Empty (Lbs)	38.2
Full Weight (Lbs)	104.7
CAPACITY	
Number of Canisters	6
Canister Dimensions (in)	4.13 x 11
1.2ml & 2.0 ml vials (5/cane)	1260
1.2ml & 2.0 ml vials (6/cane)	1512
.5 cc Straws (10/cane)	4500
.25 cc Straws (28/cane)	12600