

TW 5LD

LN₂ DEWAR

For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	5LD
LN ₂ Capacity (L)	5
Static Hold Time (Days)	6
Evap. Rate (Liters/Day)	0.77
DIMENSIONS	
Necktube (in)	5.6
Overall Dia. (in)	7.6
Height (in)	17.5
WEIGHT	
Empty (Lbs)	6.9
Full (Lbs)	15.8

FEATURES

- 5 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA