

# TW CX 100

## VAPOR SHIPPER



For the safe transport of biological samples in a vapor phase.



### FEATURES

- 3.8 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	CX100
LN <sub>2</sub> Capacity (L)	3.8
Static Hold Time (Days)	18
Evap. Rate (Liters/Day)	0.2
DIMENSIONS	
Necktube (in)	2.75
Overall Dia. (in)	8.75
Height (in)	18.5
WEIGHT	
Empty (Lbs/Kg)	11.7
Charged Weight (Lbs)	19
CAPACITY	
No. of 1.2 ml & 2.0 ml vials (5/cane)	85
No. of 1.2 ml & 2.0 ml vials (6/cane)	102
Number of .25 cc Straws (28/cane)	784
No. of .25 cc Straws, Bulk (2 levels) in Goblets	1640
Number of 1/2 straws (10/cane)	280
No. of .5 cc Straws, Bulk (2 levels) in Goblets	800
Number of Canisters	1



855.750.8191 | sales@ICBiomedical.com  
 www.ICBiomedical.com  
 Made In USA