

# TW 4D

## DRY SHIPPER

For the safe transport of biological samples in a vapor phase.



## FEATURES

- 4.1 Liter absorbed LN2 capacity
- Advanced Concept Hydrophobic Absorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	TW 4D
LN2 Capacity (L)	4.1
Static Hold Time (Days)	27
Evap. Rate (Liters/Day)	0.15
DIMENSIONS	
Necktube (in)	2.00
Overall Dia. (in)	8.75
Height (in)	18.5
CAPACITY	
No. of 1.2 ml & 2.0 ml vials (5/cane)	40
No. of 1.2 ml & 2.0 ml vials (6/cane)	48
Number of .25 cc Straws (28/cane)	448
No. of .25 cc Straws, Bulk (2 levels) in Goblets	952
Number of 1/2 straws (10/cane)	160
No. of .5 cc Straws, Bulk (2 levels) in Goblets	460
Number of Canisters	1

Note: Evaporation rate and static holding time are nominal. Actual rate may be affected by the nature of the contents, atmospheric conditions, container history and manufacturing tolerances.



855.750.8191 | [sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)  
[www.ICBiomedical.com](http://www.ICBiomedical.com)  
Made In USA