

# TW LS4800

## LABORATORY SYSTEMS



For liquid or vapor phase vial storage in convenient box-type racks.



| SPECIFICATIONS               | LS4800 |
|------------------------------|--------|
| LN <sub>2</sub> Capacity (L) | 130    |
| Static Hold Time (Days)      | 153    |
| Evap. Rate (Liters/Day)      | 0.85   |
| DIMENSIONS                   |        |
| Necktube (in)                | 8.5    |
| Overall Dia. (in)            | 26.9   |
| Height (in)                  | 35.1   |
| WEIGHT                       |        |
| Empty (Lbs)                  | 90     |
| Full (Lbs)                   | 312.4  |

## FEATURES

- Convenient Storage Systems – rack index location ring and internal location fixture, computer compatible box storage for simple inventory management
  - Number of Racks: 6
  - Number of Shelves / Rack: 8
  - 2ml Vial Capacity: 4800
- In vapor-phase storage, temperature at the top of racks is -190°C.
- High strength aluminum body and magneformed necktube design
- Superior vacuum performance with super insulation provides maximum holding times
- Security – Accessory low-level alarm and CS100 controller are available with remote connection capability
- Lockable lids



855.750.8191 | [sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)  
[www.ICBiomedical.com](http://www.ICBiomedical.com)  
Made In USA