

TW LS6000

LABORATORY SYSTEMS



For liquid- or vapor-stage vial storage in convenient box-type racks.



SPECIFICATIONS	LS6000
LN ₂ Capacity (L)	165
Static Hold Time (Days)	194
Evap. Rate (Liters/Day)	0.84
DIMENSIONS	
Necktube (in)	8.5
Overall Dia. (in)	26.9
Height (in)	39
WEIGHT	
Empty (Lbs)	121
Full (Lbs)	410

FEATURES

- Convenient Storage Systems – rack index location ring and internal location fixture, computer compatible box storage for simple inventory management
 - Number of Racks: 6
 - Number of Shelves / Rack: 10
 - 2ml Vial Capacity: 6000
- In vapor-phase storage, temperature at the top of racks is -190°C.
- High strength aluminum body and magneformed necktube design
- Superior vacuum performance with super insulation provides maximum holding times
- Security – Accessory low-level alarm and CS100 controller are available with remote connection capability
- Lockable lids



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA