

TW 12HCL

LN₂ HIGH-CAPACITY SMALL FREEZER



For long-term storage of large quantities at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	12HCL
LN ₂ Capacity (L)	12
Static Hold Time (Days)	60
Evap. Rate (Liters/Day)	0.2
DIMENSIONS	
Necktube (in/mm)	3.6 / 9.3
Overall Dia. (in/mm)	15.6/39.6
Height (in/mm)	19 / 48.3
WEIGHT	
Empty (lb/kg)	21.6 / 9.8
Full Weight (lbs/kg)	43 / 19.5
CAPACITY	
Number of Canisters	6
Canister Dimensions (in) (mm)	2.75 x 11 7 x 28
Number of .5 cc Straws – Bulk (1 Level)	2940
Number of .25 cc Straws – Bulk (1 Level)	6060
Number of 1.2 & 2.0 ml Vials	N/A



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
 Made In USA