

TW 20HC

LN₂ HIGH-CAPACITY SMALL FREEZER



For long-term storage of large quantities at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	20HC
LN ₂ Capacity (L)	20
Static Hold Time (Days)	87
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in/mm)	3.6 / 9.3
Overall Dia. (in/mm)	15.6/39.6
Height (in/mm)	24.25 / 61.6
WEIGHT	
Empty (lb/kg)	26.4 / 12
Full Weight (lb/kg)	62 / 28.1
CAPACITY	
Number of Canisters	6
Canister Dimensions (in) (mm)	2.75 x 11 7 x 28
1.2ml & 2.0 ml vials (5/cane)	570
1.2ml & 2.0 ml vials (6/cane)	684
.5 cc Straws (10/cane)	1920
.5 cc Straws- Bulk (2 Levels)	4800
.25 cc Straws (28/cane)	5376
.25 cc Straws- Bulk (2 Levels) in Goblets	10740



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
 Made In USA