

TW 35HC

LN₂ HIGH-CAPACITY SMALL FREEZER

For long-term storage of large quantities
at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid



SPECIFICATIONS	35HC
LN ₂ Capacity (L)	35
Static Hold Time (Days)	130
Evap. Rate (Liters/Day)	0.27
DIMENSIONS	
Necktube (in/mm)	4.7 / 11.9
Overall Dia. (in/mm)	18.8 / 47.8
Height (in/mm)	26.8 / 68
WEIGHT	
Empty (lb/kg)	39.0 / 17.7
Full Weight (lb/kg)	101.3 / 45.9
CAPACITY	
Number of Canisters	10
Canister Dimensions (in) (mm)	2.64 x 11 6.7 x 28
1.2ml & 2.0 ml vials (5/cane)	850
1.2ml & 2.0 ml vials (6/cane)	1020
.5 cc Straws (10/cane)	2800
.5 cc Straws- Bulk (2 Levels)	7300
.25 cc Straws (28/cane)	7280
.25 cc Straws- Bulk (2 Levels) in Goblets	16400



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA