

TW 35VHC

LN₂ HIGH-CAPACITY SMALL FREEZER

For long-term storage of large quantities
at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid



SPECIFICATIONS	35VHC
LN ₂ Capacity (L)	35
Static Hold Time (Days)	130
Evap. Rate (Liters/Day)	0.27
DIMENSIONS	
Necktube (in/mm)	4.7 / 11.9
Overall Dia. (in/mm)	18.8 / 47.8
Height (in/mm)	26.8 / 68
WEIGHT	
Empty (lb/kg)	37.9 / 17.2
Full Weight (lb/kg)	100.2 / 45.5
CAPACITY	
Number of Canisters	6
Canister Dimensions (in) (mm)	3.7 x 11 9.4 x 28
1.2ml & 2.0 ml vials (5/cane)	1050
1.2ml & 2.0 ml vials (6/cane)	1260
.5 cc Straws (10/cane)	3300
.5 cc Straws- Bulk (2 Levels)	9600
.25 cc Straws (28/cane)	9240
.25 cc Straws- Bulk (2 Levels) in Goblets	21564



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA