

# TW 38HC

## LN<sub>2</sub> STORAGE SMALL FREEZER



Designed for storing large quantities of AI/IVF samples at extremely low temperatures



### FEATURES

- Wide neck openings and large capacity provides easy access to stored samples
- Provides sample security and simplified inventory control
- Smaller footprint and lighter full shipping weight
- Numbered index rings and tapered internal canister location for ease of use
- Rugged construction - ribbed high-strength aluminum body, magenformed necktube design and durable paint
- 5-Year vacuum warranty

SPECIFICATIONS	38HC
LN <sub>2</sub> Capacity (L)	38
Static Hold Time (Days) Liquid	120
Working Time Days (Liquid)	75
Evap. Rate (Liters/Day)	0.31
DIMENSIONS	
Neck Diameter (in/mm)	5.0 / 127
Overall Height (in/mm)	26.8 / 680.7
Overall Diameter (in/mm)	18.8 / 477.5
WEIGHT	
Weight Empty (lb/kg)	38.2 / 17.3
Weight Full (lb/kg)	105.9 / 48
CAPACITY	
Number of Canisters	10
Canister Dimensions (in) (mm)	2.75 x 11 69.9 x 279.4
# of 1.2 & 2.0 ml Vials (5/Cane)	950
# of 1.2 & 2.0 ml Vials (6/Cane)	1140
# of .05 cc Straws (10/Cane)	3200
# of .05 cc Straws Bulk (2 Levels)	8000
# of .25 cc Straws (28/Cane)	8960
# of .25 cc Straws Bulk (2 Levels) in Goblet	17,200



855.750.8191 | [sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)

[www.ICBiomedical.com](http://www.ICBiomedical.com)

Made In USA