

TW 38VHC

LN₂ HIGH-CAPACITY SMALL FREEZER

For long-term storage of large quantities
at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid



SPECIFICATIONS	38VHC
LN ₂ Capacity (L)	38
Static Hold Time (Days)	120
Evap. Rate (Liters/Day)	0.31
DIMENSIONS	
Necktube (in/mm)	5.00 / 127
Overall Dia. (in mm)	18.8 / 477.5
Height (in/mm)	26.8 / 680.7
WEIGHT	
Empty (lb/kg)	38.2 / 17.3
Full Weight (lb/kg)	104.7 / 47.5
CAPACITY	
Number of Canisters	6
Canister Dimensions (in) (mm)	4.13 x 11 104.9 x 279.4
1.2ml & 2.0 ml vials (5/cane)	1260
1.2ml & 2.0 ml vials (6/cane)	1512
.5 cc Straws (10/cane)	4500
.25 cc Straws (28/cane)	12600



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA