

**TW 38VHC**
**LN<sub>2</sub> HIGH-CAPACITY  
SMALL FREEZER**

For long-term storage of large quantities at cryogenic temperatures.



## FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	38VHC
LN <sub>2</sub> Capacity (L)	38
Static Hold Time (Days)	120
Evap. Rate (Liters/Day)	0.31
DIMENSIONS	
Necktube (in/mm)	5.00 / 12.7
Overall Dia. (in mm)	18.8 / 47.8
Height (in/mm)	26.8 / 68
WEIGHT	
Empty (lb/kg)	38.2 / 17.3
Full Weight (lb/kg)	104.7 / 47.5
CAPACITY	
Number of Canisters	6
Canister Dimensions (in) (mm)	4.13 x 11 10.5 x 28
1.2ml & 2.0 ml vials (5/cane)	1260
1.2ml & 2.0 ml vials (6/cane)	1512
.5 cc Straws (10/cane)	4500
.25 cc Straws (28/cane)	12600