

10LD

LN₂ DEWAR

For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	10LD
LN ₂ Capacity (L)	10
Static Hold Time (Days)	45
Evap. Rate (Liters/Day)	0.22
DIMENSIONS	
Necktube (in/mm)	2.0 / 50.8
Overall Dia. (in/mm)	11.4 / 289.6
Height (in/mm)	23.5 / 597
WEIGHT	
Empty (lb/kg)	14.5 / 6.6
Full (lb/kg)	32.3 / 14.7

FEATURES

- 10 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA