

25LD

LN₂ DEWAR

For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	25LD
LN ₂ Capacity (L)	25
Static Hold Time (Days)	109
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in/mm)	2.5 / 63.5
Overall Dia. (in/mm)	15.6 / 396.2
Height (in/mm)	25.8 / 655.3
WEIGHT	
Empty (lb/kg)	23.2 / 10.5
Full (lb/kg)	67.7 / 30.7

FEATURES

- 25 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency