

12HCL

LN₂ HIGH-CAPACITY SMALL FREEZER

For long-term storage of large quantities
at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid



SPECIFICATIONS	12HCL
LN ₂ Capacity (L)	12
Static Hold Time (Days)	60
Evap. Rate (Liters/Day)	0.2
DIMENSIONS	
Necktube (in/mm)	3.6 / 91.4
Overall Dia. (in/mm)	15.6/396.2
Height (in/mm)	19 / 482.6
WEIGHT	
Empty (lb/kg)	21.6 / 9.8
Full Weight (lbs/kg)	43 / 19.5
CAPACITY	
Number of Canisters	6
Canister Dimensions (in) (mm)	2.75 x 11 70 x 279.4
Number of .5 cc Straws – Bulk (1 Level)	2940
Number of .25 cc Straws – Bulk (1 Level)	6060
Number of 1.2 & 2.0 ml Vials	N/A



855.750.8191 | sales@ICBiomedical.com

www.ICBiomedical.com

Made In USA

© 2022, IC Biomedical, LLC