

34HC

LN₂ HIGH-CAPACITY SMALL FREEZER

For long-term storage of large quantities at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	34HC
LN ₂ Capacity (L)	34
Static Hold Time (Days)	200
Evap. Rate (Liters/Day)	0.17
DIMENSIONS	
Necktube (in/mm)	3.6 / 91.4
Overall Dia. (in/mm)	18.8 / 477.5
Height (in/mm)	26.3 / 66.8
WEIGHT	
Empty (lb/kg)	35.38 / 16
Full Weight (lb/kg)	95.9 / 43.5
CAPACITY	
Number of Canisters	6
Canister Dimensions (in) (mm)	2.75 x 11 70 x 279.4
1.2ml & 2.0 ml vials (5/cane)	570
1.2ml & 2.0 ml vials (6/cane)	684
.5 cc Straws (10/cane)	1920
.5 cc Straws- Bulk (2 Levels)	4800
.25 cc Straws (28/cane)	5376
.25 cc Straws- Bulk (2 Levels) in Goblets	10740