

# 4LD

## LN<sub>2</sub> DEWAR

For safe storage and dispensing of small amounts of liquid nitrogen.



| SPECIFICATIONS               | 4LD        |
|------------------------------|------------|
| LN <sub>2</sub> Capacity (L) | 4          |
| Static Hold Time (Days)      | 10         |
| Evap. Rate (Liters/Day)      | 0.40       |
| DIMENSIONS                   |            |
| Necktube (in/mm)             | 1.2 / 30.5 |
| Outside Dia. (in/mm)         | 7.6 / 193  |
| Height (in/mm)               | 17 / 432   |
| WEIGHT                       |            |
| Empty (lb/kg)                | 6.6 / 3.0  |
| Full (lb/kg)                 | 13.7 / 6.2 |

## FEATURES

- 4 Liter LN<sub>2</sub> capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency